

Application Memo

Concentration of Acetic Acid

Industry	Food
Instrument	Automatic potentiometric titrator
Measurement method	Acid-base titration
Standards	

1. Overview

In this memo, the concentration of the prepared 2w/w% acetic acid sample is tested. The diluted sample is titrated with 0.1mol/L sodium hydroxide up to pH8.3. The acetic acid concentration is calculated from the titration volume of sodium hydroxide.


2. Apparatus

Main unit	Automatic potentiometric titrator (preamplifier STD)
Electrode	Combined glass electrode Temperature compensation electrode

3. Reagents

Titrant	0.1mol/L sodium hydroxide
Additive	Pure water

4. Example



—Titration curve—

—Measurement results—			
	Sample (g)	Titer (mL)	Concentration (w/w%)
1	1.492	5.0831	1.9787
2	1.492	5.0745	1.9754
3	1.492	5.0814	1.9781
Average			1.9774
SD			0.0018
RSD(%)			0.089

Please feel free to contact us for any further information.
 <Contact> Kyoto Electronics Manufacturing Co., Ltd.
 Overseas Sales & Marketing Sect.
<http://www.kyoto-kem.com/en/contact/form.php>