

Application Note

Measuring the acidity of fruit beverages using a multiple sample changer

Industry	Food & beverage
Instrument	Automatic potentiometric titrator
Measurement method	Potentiometric titration / Neutralization titration
Standards	-

1. Scope

Fruit beverages contain citric acid, malic acid, and other organic acids. The acidity indicates the amount of these organic acids. In fruit beverages, the ratio of acidity to sugar content has a significant effect on the preference and flavor. Accordingly, it is very important to control the acidity in order to manufacture beverages of consistent quality.

This Application Note describes an example of the measurement of the acidity of a fruit beverage (from mandarin oranges). Note that the measurements were performed using a multiple sample changer. Multiple beakers containing the samples were placed on a turntable, enabling consecutive measurements to be performed automatically. In this test, the measurements were performed with the parameter settings configured so that a pH of 8.0 was detected as the end point.

2. Post-measurement procedure

To suppress efflux and enrichment of the electrolyte when storing the electrodes, seal the electrolyte filling port in the combined glass electrode with a rubber stopper.

The performance of the combined glass electrode quickly deteriorates if it is stored while dry. The following storage methods are recommended.

- Short-term storage (less than one month): Store in pure water.
- Long-term storage (one month or longer): Store in a mixture of pH4 standard solution and 3.3 mol/L potassium chloride solution in a 1:1 volume ratio.

3. Apparatus

- Automatic potentiometric titrator (preamplifier STD)
- Multiple sample changer
- Combined glass electrode (Electrolyte: 3.3 mol/L potassium chloride solution)
- Temperature compensated electrode

4. Reagents

Titration liquid 0.1 mol/L sodium hydroxide solution

Others Standard solutions at a pH of 7 and 9

5. Procedure

- 1) Calibrate the combined glass electrode using the standard solutions at pH 7 and 9.
- 2) Introduce the sample into a beaker and weigh it.
- 3) Using a measuring cylinder, add 100 mL of pure water, and place the beaker in the multiple sample changer.
- 4) Titrate with a 0.1 mol/L sodium hydroxide solution, with a pH of 8.0 as the end point.

6. Calculation

$$\text{Acidity (\%)*} = (\text{EP1} - \text{BL1}) \times \text{TF} \times 6.404 \times 0.1/\text{S}$$

EP1 Titration volume (mL)

BL1 Titration volume (mL) for a blank test = 0

TF Titration liquid factor = 0.9920

S Amount of sample introduced (g)

* Calculated as the mass percent concentration of citric acid

7. Example

—Titration parameter—

<u><Titr. Mode></u>	Auto Int.	<u><Ctrl. Para.></u>	
<u><Titr. Form></u>	Level Stop	Number of EP	1
		1st End Level	8.0 pH
<u><Titr. Para.></u>		Gain	1
Max Volume	20(mL)	Data Sampling	Auto
Channel/Unit(Ctrl.)	Ch1, pH	Ctrl. Speed	Standard
Channel/Unit(Ref.)	Off	Other Control	Standard
pH Polarity	Standard	Auto Int. Mode	Standard
Titr. Type Check	No Check	Stirrer Speed	4
Direction	Auto		
Wait Time	0 (s)		
Dose Mode	None		
<u><Sample Changer></u>			
Shower	15s		
Drain	30s		

(The above condition is an example. The setting condition depends on the model.)

—Titration curve—




Table 1 Measurement result

	Sample (g)	Titration amount (mL)	Acidity (%)
1	10.1474	9.9927	0.626
2	10.0642	9.9089	0.625
3	10.2621	10.0971	0.625
4	10.1302	9.9625	0.625
5	10.5239	10.3445	0.624
6	10.1199	9.9523	0.625
7	10.4572	10.2788	0.624
8	10.2777	10.1057	0.625
9	10.7381	10.5577	0.625
10	10.1543	9.9836	0.625
Mean	-	-	0.625
SD	-	-	0.0004
RSD(%)	-	-	0.065

8. Summary

Acidity measurements were performed using an automatic potentiometric titrator. In this test, a multiple sample changer was used to perform a consecutive series of measurements automatically. Using a multiple sample changer makes it possible to reduce the labor required for the measurements and improve efficiency.