


Application Note

Ion exchange capacity of cation exchange resin

Industry	Chemicals
Instrument	Automatic potentiometric titrator
Measurement method	Potentiometric titration / Neutralization titration
Standards	

1. Scope

This Application Note describes an example of measuring the ion exchange capacity of a strongly acidic hydrogen ion type cation exchange resin. A strongly acidic hydrogen ion type cation exchange resin contains the sulfo group (-SO₃H), and as indicated below, exchanges sodium ions, calcium ions, and other cations in the solvent with hydrogen ions.



In terms of measuring the ion exchange capacity, in light of the properties mentioned above, an applicable method is to add a sodium chloride solution to the cation exchange resin, and then titrate the hydrogen ions released with a sodium hydroxide solution. The exchange capacity is calculated as (mmol/g).

2. Post-measurement procedure

To suppress efflux and enrichment of the electrolyte when storing the electrodes, seal the electrolyte filling port in the combined glass electrode with a rubber stopper.

The performance of the combined glass electrode quickly deteriorates if it is stored while dry. The following storage methods are recommended.

- Short term storage (less than one month): Store it submerged in pure water.
- Long term storage (at least one month): Store it submerged in a 1:1 volume ratio mixture of a pH 4 standard liquid and a 3.3 mol/L aqueous potassium chloride solution.

3. Apparatus

Main unit	Automatic potentiometric titrator (standard preamplifier STD)
Electrode	Combined glass electrode (Electrolyte 3.3 mol/L aqueous potassium chloride solution) Temperature compensated electrode

4. Reagents

Titration liquid	0.1 mol/L sodium hydroxide solution
Additive reagent	1 mol/L sodium chloride solution

5. Procedure

- 1) Introduce 0.5 g of the cation exchange resin into a beaker and weigh it.
- 2) Add 10 mL of a 1 mol/L aqueous sodium chloride solution and stir for 1 minute.
- 3) Add pure water to dilute the solution up to approximately 100 mL, and titrate with a 0.1 mol/L aqueous sodium hydroxide solution.

6. Calculation

$$\text{Exchange capacity (mmol/g)} = (\text{EP1} - \text{BL1}) \times \text{TF} \times \text{C1} \times \text{K1/S}$$

EP1	Titration volume required to the first endpoint (mL)
BL1	Titration volume required for a blank test (0 mL)
TF	Titration solution factor (1.0003)
C1	Concentration conversion coefficient (0.0001 mol/mL)
K1	Unit conversion coefficient (1000)
S	Amount of sample introduced (g)


7. Example

— Parameter —

<u><Titr. Mode></u>	Auto Int.	<u><Ctrl. Para.></u>	
<u><Titr. Form></u>	Level Stop	Number of EP	1
		End Sense	Auto
<u><Titr. Para.></u>		Gain	1
Max Volume	20 (mL)	Data Sampling	Auto
Channel/Unit(Ctrl.)	Ch1, mV	Ctrl. Speed	Standard
Wait Time	0 (s)	Other Control	Standard
Dose Mode	None	Stirrer Speed	4

(The above condition is an example. The setting condition depends on the model.)

— Example of Titration curve —



— Measurement results —

Table 1 Measurement result of exchange capacity of cation exchange resin

	Sample (g)	Titration (mL)	Ion exchange capacity (mmol/g)
1	0.4940	6.2339	1.2623
2	0.4859	6.1134	1.2585
3	0.4911	6.2120	1.2653
Mean			1.2620
SD			0.0034
RSD (%)			0.27