

Application Note

Determination of inorganic sulfate in bioethanol

Industry	Chemicals
Instrument	Automatic potentiometric titrator
Measurement method	Precipitation titration
Standards	ASTM D7318

1. Scope

Bioethanol for fuel is sometimes contaminated with inorganic sulfates, which can cause engine deterioration and air pollution. Therefore, it is necessary to confirm that the inorganic sulfate is less than a specified amount. This Application Note introduces an example of measuring sulfate ion concentration in ethanol using the method specified by ASTM D7318.

The quantification method for this standard is shown below.

In a solution of 100 g of sulfate-free ethanol, add 0.20 mL of 0.01 mol/L aqueous sulfate standard solution, and use the titration volume of this solution as a blank test value. Add the same volume of the aqueous sulfate standard to 100 g of ethanol sample. The difference in titration volume from the blank test is defined as the concentration of sulfate in the sample.

2. Precautions

- 1) Before starting measurements each day, purge the titrant several times between the reagent bottle and the burette to equalize the titrant concentration, then discharge about 10 mL of the titrant between the burette and the titration nozzle.
- 2) This measurement uses the reagents harmful to the human body. When you follow this application note, wear masks, gloves, protective equipment, etc.

3. Apparatus

Main unit	Automatic potentiometric titrator (Preamplifier STD)
Electrode	Lead ion electrode, Double junction reference electrode
Outer cylinder inner liquid	1 mol/L lithium chloride solution (ethanol solvent)

4. Reagents

Titrant	0.025 mol/L Lead nitrate standard solution
Additive reagents	0.01 mol/L Aqueous sulfate standard
	0.1 mol/L Perchloric acid solution
	Ethanol* (for blank test) *Used a sulfate-free reagent

5. Procedure

-Blank test-

- 1) Collect 100 g of ethanol (sulfate-free) and measure the mass.
- 2) Add 0.20 mL of 0.01 mol/L aqueous sulfate standard.
- 3) Add 1 mL of 0.1 mol/L perchloric acid solution. Check the pH with pH test paper (range of 3 to 5). If the sample exceeds pH 5, add 0.1 mol/L perchloric acid solution as necessary, and adjust to a pH of 3 to 5.
- 4) Titrate with 0.0025 mol/L lead nitrate standard solution.

-Measurement-

- 1) Collect 100 g of ethanol sample and measure the mass.
- 2) Add 0.20 mL of 0.01 mol/L aqueous sulfate standard.
- 3) Add 1 mL of 0.1 mol/L perchloric acid solution. Check the pH with pH test paper (range of 3 to 5). If the sample exceeds pH 5, add 0.1 mol/L perchloric acid solution as necessary, and adjust to a pH of 3 to 5.
- 4) Titrate with 0.0025 mol/L lead nitrate standard solution.

6. Calculation

$$\text{Sulfate ion concentration (mg/kg)} = (\text{EP1} - \text{BL1}) \times \text{TF} \times \text{C1} \times \text{K1/S}$$

EP1	Titration amount (mL)
BL1	Titration amount (mL) of blank test = 0.7955 (mL)
TF	Factor of titrant = 1.0259
C1	Concentration conversion coefficient = 0.24014
K1	Unit conversion factor = 1000
S	Sample size (g)

7. Example

— Parameter —

<Titration Mode>	Auto Intermit	<Control Parameter>	
		Number of EP	1
<Titration Form>	EP Stop	End Sense	Set
		Gain	3
<Titration Parameter>		Data Sampling	Standard
Channel, Unit(For Control)	Ch1, mV	Control Speed Mode	Standard
Channel, Unit(For Reference)	Off	Other Control	Standard
pH Polarity	Standard	Auto Int. Mode	Standard (Sample) Blank (Blank)
Type of Titration	Not Check.	Stirrer Speed	4
EP Direction	Auto		
Wait Time	30 (s)		
Dose Mode	Off		

(The above condition is an example. The setting condition depends on the model.)


— Measurement results —

The measurement results for each sample are shown in Tables 1 and 2.

	Sample (g)	Titration (mL)	Conc. (mg/kg)
1	100.0141	1.0527	0.63
2	100.0755	1.0303	0.58
3	100.0451	1.0490	0.62
Average	-	-	0.61
SD	-	-	0.03
RSD (%)	-	-	4.3

	Sample (g)	Titration (mL)	Conc. (mg/kg)
1	100.0603	1.2495	1.12
2	100.0008	1.2220	1.05
3	100.0255	1.2469	1.11
Average	-	-	1.09
SD	-	-	0.04
RSD (%)	-	-	3.5

— Example of titration curve —



Example of titration curve

8. Reference

ASTM D7318-19 Standard Test Method for Existent Inorganic Sulfate in Ethanol by Potentiometric Titration