


Application Note

Standardization of potassium bromide-potassium bromate standard solution

| | |
|--------------------|---------------------------------------|
| Industry | Petrochemicals |
| Instrument | Automatic potentiometric titrator |
| Measurement method | Redox titration |
| Standards | UOP Method 304, ASTM D1159, JIS K2605 |

1. Scope

Potassium bromide-potassium bromate standard solution is used as a titrant for measuring bromine number and bromine index. This application note describes the standardization of 0.25mol/L potassium bromide potassium bromate standard solution according to the UOP Method 304. These reactions involved in the standardization are outlined below.



2. Precautions

- Prior to the first titration, prepare the burette by dispensing and aspirating the titrant between it and the reagent bottle several times. Once done, dispense approximately 10mL of the titrant.
- The inner electrolyte of the combined platinum electrode should be changed every two weeks.

3. Apparatus

| | |
|-----------|------------------------------------------------------------------------------------------|
| Main unit | Automatic potentiometric titrator (Preamplifier STD) |
| Electrode | Combined platinum electrode (Inner electrolyte: 3.3mol/L Potassium chloride solution) |

4. Reagents

| | |
|-------------------|-----------------------------------------------------|
| Titrant | 0.1mol/L Sodium thiosulfate solution |
| Additive reagents | 15% Potassium iodide solution |
| Reagents | Glacial acetic acid, Concentrated hydrochloric acid |

5. Procedures

- 1) Combine 50 mL of glacial acetic acid and 1 mL of concentrated hydrochloric acid into a glass stoppered Erlenmeyer flask, and cool in an ice bath for 10 minutes while stirring.
- 2) Add 5 mL of 0.25 mol/L potassium bromide-potassium bromate solution, 1 or 2 drops per second while stirring.
- 3) Stopper the flask and cool in an ice bath for 5 minutes.
- 4) Add 5mL of 15% KI solution and shake vigorously.
- 5) Transfer the reaction mixture to a beaker, rinse the inside of the flask and glass stopper with 100mL of pure water and then transfer this rinse water to the same beaker.
- 6) Titrate the solution with 0.1mol/L sodium thiosulfate solution.

6. Calculation

$$\text{Factor} = \text{EP1} \times \text{C1} \times \text{FA} / \text{S}$$

| | |
|-----|----------------------------------------------------------|
| EP1 | Titration amount (mL) |
| C1 | Concentration conversion coefficient (0.2) |
| FA | Factor of 0.1 mol/L sodium thiosulfate solution (1.0043) |
| S | Quantity of sample (mL) |

7. Example

-Parameter-

| <u><Titr. Mode></u> | Auto Int. | <u><Ctrl. Para.></u> | |
|----------------------------|-------------|----------------------------|----------|
| <u><Titr. Form></u> | EP Stop | Number of EP | 1 |
| <u><Titr. Para.></u> | | End Sense | Auto |
| Max. Volume (mL) | 30 | Gain | 1 |
| Channel / Unit (Ctrl.) | Ch1, mV | Data Sampling | Auto |
| Channel / Unit (Ref.) | Off | Ctrl. Speed Mode | Standard |
| pH Polarity | Standard | Other Control | Standard |
| Type of Titration | Not check | Auto Int. Mode | Standard |
| EP Direction | Auto | Stirrer Speed | 4 |
| Wait Time (s) | 0 | | |
| Dose Mode | Volume Stop | | |
| Stop Volume (mL) | 20 | | |
| Cut Off Time | 0 | | |
| Disp. Speed (s/mL) | 1 | | |
| Wait Time (s) | 0 | | |

(Listed above are example settings. Availability of settings may vary by instrument model.)

-Example of titration curve-




Table 1 Standardization of 0.25mol/L potassium bromide - potassium bromate solution

| | Titration (mL) | Factor |
|---------|----------------|--------|
| 1 | 25.1950 | 1.0121 |
| 2 | 25.1841 | 1.0117 |
| 3 | 25.1726 | 1.0112 |
| Mean | - | 1.0117 |
| SD | - | 0.0005 |
| RSD (%) | - | 0.05 |

8. Reference

UOP Method 304 BROMINE NUMBER AND BROMINE INDEX OF HYDROCARBONS BY POTENTIOMETRIC TITRATION