

Application Note

Ferrous ion and ferric ion in pickling solution

Industry	:	Iron and steel
Instrument	:	Automatic potentiometric titrator
Measurement method	:	Photometric titration / Chelatometric titration
Standards	:	

1. Scope

Ferrous ion (Fe^{2+}) and ferric ion (Fe^{3+}) in pickling solution were measured by chelatometric titration. **Combined method*** enables measurement of Fe^{2+} and Fe^{3+} with one sample.

First, a sample was titrated with 0.1 mol/L ethylenediaminetetraacetic acid · 2Na (EDTA), using salicylic acid as an indicator. An inflection point on the titration curve was regarded as the end point, and Fe^{3+} concentration was calculated from the volume of EDTA solution consumed to titrate sample to the end point. Then, Fe^{2+} in the solution after titration was oxidized to Fe^{3+} by adding potassium peroxodisulfate, and the solution was titrated with 0.1 mol/L EDTA solution. Total concentration of Fe^{2+} and Fe^{3+} was calculated from the volume of EDTA solution consumed to the end point. Fe^{2+} concentration was calculated by subtracting the Fe^{3+} concentration from the total concentration of Fe^{2+} and Fe^{3+} .

*Combined method

Combined method is used for perform multiple measurement on one sample. It is possible to connect up to 5 methods.

2. Precautions

- 1) Adjust pH of a titration solution to 2-3 before measurement.
- 2) Handle the reagents in a well ventilated room or a draft chamber.

3. Post-measurement procedure

Wash photometric sensor with ethanol and then water.

4. Apparatus

Main unit : Automatic potentiometric titrator (preamplifier : PTA)
Electrode : Photometric sensor (interference filter : 530 nm)

5. Reagents

Titrant : 0.1 mol/L EDTA solution
Indicator : 2% salicylic acid ethanol solution
Addition reagent : potassium peroxodisulfate

6. Procedure

-Measurement-

- Method 1 ↑
- Method 2 ↓
- 1) Add 5 mL of 10-fold diluted sample into a 200 mL tall beaker.
 - 2) Add 100 mL of pure water and 1 mL of 2% salicylic acid ethanol solution.
 - 3) Titrate with 0.1 mol/L EDTA solution to measure Fe³⁺ concentration.
 - 4) Add 1 g of potassium peroxodisulfate into the beaker and dissolve it by stirrer during “Wait time” of combined method 2. After dissolving, press [Start] button to skip “Wait time” and start titration.
 - 5) Titrate with 0.1 mol/L EDTA solution to measure total concentration of Fe²⁺ and Fe³⁺.
 - 6) Calculate Fe²⁺ concentration by subtracting the Fe³⁺ concentration from the total concentration of Fe²⁺ and Fe³⁺.

7. Calculation

Combined method 1 : Measurement of Fe³⁺

$$\text{Fe}^{3+} \text{ (g/L)} = (\text{EP1} - \text{BL1}) \times \text{TF} \times \text{C1} \times \text{K1} / (\text{R} \times \text{S})$$

EP1	:	Titer (mL)
BL1	:	Blank level = 0.0000 mL
TF	:	Factor of titrant = 0.9993
C1	:	Concentration conversion coefficient = 5.585 mg/mL
K1	:	Unit conversion coefficient = 1
R	:	Dilution factor = 0.1
S	:	Quantity of diluted sample (mL)

Combined method 2 : Measurement of Fe²⁺

$$\text{Fe}^{2+} \text{ and Fe}^{3+} \text{ (g/L)} = (\text{EP1} - \text{BL1}) \times \text{TF} \times \text{C1} \times \text{K1} / (\text{R} \times \text{S})$$

EP1	:	Titer (mL)
BL1	:	Blank level = 0.0000 mL
TF	:	Factor of titrant = 0.9993
C1	:	Concentration conversion coefficient = 5.585 mg/mL
K1	:	Unit conversion coefficient = 1
R	:	Dilution factor = 0.1
S	:	Quantity of diluted sample (mL)

$$\text{Fe}^{2+} \text{ (g/L)} = \text{CO1} - \text{FCO1}$$

CO1	:	Total concentration of Fe ²⁺ and Fe ³⁺ (g/L)
FCO1	:	Fe ³⁺ (g/L)

8. Example

-Titration parameter-

Combined method

<Titr. Mode> : Auto Intermit

<Titr. Form> : EP Stop

<Titr. Para.>

Max. Volume : 20 (mL)

Channel/Unit(Ctrl.) : Ch3, %T

Channel/Unit(Ref.) : Off

pH Polarity : Standard

Tit. Type Check : No Check

Direction : Auto

Wait Time : 0 (s) (Combined method 1)

: 600 (s) (Combined method 2)

Dose Mode : None

<Ctrl. Para.>

Number of EP : 1

End Sense : Auto

Gain : 1

Data Sampling : Auto

Ctrl. Speed : Standard


Other Ctrl. : Standard

Auto Int. Mode : Standard

Stirrer Speed : 4

(The measurement parameter and the titration curve are an example of our automatic potentiometric titrator. For other models, parameter item may be different or other parameter item may be added.)

-Results-



-Measurement results-

	Quantity of diluted sample (mL)	Method 1 titer (mL)	Fe ³⁺ (g/L)	Method 2 titer (mL)	Fe ²⁺ and Fe ³⁺ (g/L)	Fe ²⁺ (g/L)
1	5	1.7639	19.69	4.6807	52.25	32.56
2	5	1.7823	19.89	4.6907	52.36	32.46
3	5	1.7650	19.70	4.6810	52.25	32.55
Mean	-	-	19.76	-	52.29	32.52
SD	-	-	0.11	-	0.06	0.06
RSD (%)	-	-	0.57	-	0.12	0.17

9. Summary

In this measurement, the results showed a good repeatability with below 1% RSD (Relative standard deviation), and continuous titration of Fe³⁺ and total iron ion with “combined method” can be performed. When quantity of Fe²⁺ is less than the quantity equivalent to excess quantity of titrant in “method 1”, continuous titration by “combined method” cannot be performed. In that case, titrate Fe³⁺ and total iron ion separately.

In some samples, verification of the measurement capability is required. In such case, please contact us.

10. References