

Application Note

Density of AdBlue®

Industry	:	Chemicals
Instrument	:	Density / specific gravity meter
Measurement method	:	Resonant frequency oscillation
Standards	:	ISO 22241, JIS K2247

1. Scope

Density of AdBlue® was determined based on “ISO 22241 Diesel engines - NOx reduction agent AUS 32 -”.

2. Precautions

- 1) Use sufficiently dried desiccant.
- 2) Perform a calibration with dry air and pure water before measurement.

3. Post-measurement procedure

- 1) Drain a sample in the measuring cell and wash it with pure water. Then rinse the cell with ethanol and dry it sufficiently.

4. Apparatus

Main unit : Density / specific gravity meter

5. Reagents

Rinse liquid (for washing) : Pure water
Rinse liquid (for drying) : Ethanol

6. Procedure

-Preparation-

- 1) Set the measuring temperature of the main unit to 20°C, and wait until it is stabilized.
- 2) Perform a calibration with dry air and pure water.

-Measurement-

- 1) Introduce a sample into the measuring cell without letting air bubbles enter into the cell.
- 2) Start a measurement.
- 3) After a measurement, wash the cell with pure water and ethanol, and dry it.

7. Example

-Parameter-

<Measurement Parameter>

Set temperature : 20.00 (°C)
 Stability sense : 1
 Limit time : 600 (sec)
 Viscosity correction : 0 (off)
 Calibration material : 0 (air and water)

<Contents>

Display item : Density (g/cm³)
 Decimal place : 4

<Temperature Compensation>

Temp. Comp. : 0 (off)

(This parameter is an example of our refractometer. For other models, parameter items may be different or other items may be added.)

-Measurement results-

	Density (g/cm ³)
1	1.0905
2	1.0905
Mean	1.0905
Repeatability*	0.0000

*Repeatability means the difference of the two results.

8. Summary

As the measurement results of density of AdBlue®, the repeatability was within the limits (ISO 3675 : 0.0005, ISO 12185 : 0.0002) specified by the standards, and the results satisfied the quality requirement (within 1.0870 - 1.0930 g/cm³).

When actually measuring, please refer to the latest standards.

9. References

- 1) ISO 22241 : 2006 (E) Diesel engines — NOx reduction agent AUS 32 —
- 2) JIS K2247 : 2009 Diesel engines — NOx reduction agent AUS 32 —