



# **Application Memo Purity of Benzalkonium Chloride**

Industry Pharmaceutical

Instrument Automatic potentiometric titrator

Measurement method Redox titration

Standards Japanese pharmacopoeia codex 17

#### 1. Overview

Benzalkonium chloride (Ammonium salt, class 4) is a bactericide. Here, its content (quantified as  $C_{22}H_{40}ClN=354.02$ ) in a sample from commercially sold bactericide is measured according to Japanese pharmacopoeia codex 17 by potentiometric titration with 0.02mol/L sodium tetraphenylboron solution. The endpoint is determined by the maximum inflexion point on the titration curve.

## 2. Apparatus

Main unit Automatic potentiometric titrator (preamplifier STD)


Detector Combined platinum electrode

### 3. Reagents

Titrant 0.02mo1/L sodium tetraphenylboron solution

Additive Diluted hydrochloric acid (1+2)

## 4. Example



-Measurement results-			
	Sample	Titer	Conc.
	(g)	(mL)	(%)
1	0.1589	20.4143	95.51
2	0.1536	19.6518	95.11
3	0.1493	19.2010	95.61
Average			95.41
SD			0.26
RSD(%)			0.27

Please feel free to contact us for any further information.

< Contact > Kyoto Electronics Manufacturing Co., Ltd.

Overseas Sales & Marketing Sect.

http://www.kyoto-kem.com/en/contact/form.php

