

Application Memo

Polyvinyl Potassium Sulfate (PVSK)

Industry	Cosmetics soap
Instrument	Automatic potentiometric titrator
Measurement method	Colloidal titration
Standards	

1. Overview

Colloidal titration is a polyelectrolyte titration method devised in 1946. Its principle is based on the immediate formation of complex by polycation and polyanion. The endpoint is usually found by the color change of the toluidine blue indicator, but this time we will use the particle changer detector unit PCD-500. This unit can automatically detect the end point without an indicator.

The test sample is titrated with 0.0005mol/L methyl glycol chitosan by PCD-500, and the concentration of polyvinyl potassium sulfate is calculated from the titration volume up to the endpoint.


2. Apparatus

Main unit Automatic potentiometric titrator (preamplifier STD)
Particle changer detector unit (PCD-500)

3. Reagents

Titrant 0.0005mol/L methyl glycol chitosan
Solvent Pure water

4. Example



—Titration curve—

—Measurement results—

	Sample (mL)	Titer (mL)	PVSK (mol/L)
1	5.000	4.7759	0.4776
2	5.000	4.8163	0.4816
3	5.000	4.7617	0.4762
Average			0.4785
SD			0.0028
RSD(%)			0.59

Please feel free to contact us for any further information.
<Contact> Kyoto Electronics Manufacturing Co., Ltd.
Overseas Sales & Marketing Sect.
<http://www.kyoto-kem.com/en/contact/form.php>