


Application Memo

Cationic Surfactant in Rinse Agent

Industry	Cosmetics soap
Instrument	Automatic potentiometric titrator
Measurement method	Ion association titration
Standards	

1. Overview

Cationic surfactant concentration can be calculated from the potentiometric titration volume using an electrode for surface-active agent. Titrate with 0.005mol/L dodecyl sodium sulfate (SDS) up to the endpoint. The endpoint is the maximum inflexion point on the titration curve.



Where : R CH₃(CH₂)_nO (n=9~13)
R₄ Carbon chain of surfactant


2. Apparatus

Main unit	Automatic potentiometric titrator (preamplifier STD)
Electrode	Combined surfactant electrode

3. Reagents

Titrant	0.005mol/L SDS
Solvent	1% TritonX-100, 0.01mol/L hydrochloric acid

4. Example



—Measurement results—

	Sample (g)	Titer (mL)	Concentration (%)
1	5.3884	3.0665	0.7887
2	5.3884	3.0712	0.7899
3	5.3884	3.1444	0.8087
Average			0.7958
SD			0.0112
RSD(%)			1.41

Please feel free to contact us for any further information.
 < Contact > Kyoto Electronics Manufacturing Co., Ltd.
 Overseas Sales & Marketing Sect.
<http://www.kyoto-kem.com/en/contact/form.php>