

Application Catalog

Acid Number of Oil in Lipstick

| | |
|--------------------|-----------------------------------|
| Industry | Cosmetics & soap |
| Instrument | Automatic potentiometric titrator |
| Measurement method | Acid-base titration |
| Standards | |

1. Overview

A test sample is dissolved in mixture of toluene, 2-propanol and a small amount of water, and the sample is titrated by potentiometry with 0.1mol/L potassium hydroxide + 2-propanol solution. The endpoint is obtained on the titration curve.

The acid number is calculated from the titration volume up to the endpoint.


2. Apparatus

| | |
|-----------|---|
| Main unit | Automatic potentiometric titrator (preamplifier STD) |
| Electrode | pH Glass electrode Double junction reference electrode Temperature compensation electrode |

3. Reagents

| | |
|---------|--|
| Titrant | 0.1mol/L potassium hydroxide + 2-propanol solution |
| Solvent | Toluene, Pure water, 2-propanol |

4. Example



-Titration curve-

-Measurement results-

| | Sample (g) | Titer (mL) | Acid number (mg/g) |
|---------|---------------|---------------|--------------------------|
| 1 | 1.0040 | 0.5904 | 2.488 |
| 2 | 1.0038 | 0.5666 | 2.355 |
| 3 | 1.0858 | 0.6606 | 2.663 |
| Average | | | 2.502 |
| SD | | | 0.154 |
| RSD(%) | | | 6.17 |

Please feel free to contact us for any further information.

< Contact > Kyoto Electronics Manufacturing Co., Ltd.

Overseas Sales & Marketing Sect.

<http://www.kyoto-kem.com/en/contact/form.php>