

## Application Memo

# Copper Ion (Cu<sup>2+</sup>) in Sulfuric Acid (20g/L)

Industry	Inorganic chemical industry
Instrument	Automatic potentiometric titrator
Measurement method	Redox titration
Standards	

## 1. Overview

Here we demonstrate quantification of copper in sulfuric acid by potentiometric redox titration. First, add pure water and 10% potassium iodide to the sample liquid, and titrate with the 0.1mol/L sodium thiosulfate. The endpoint is the max inflexion on the titration curve. The concentration of copper is calculated from the titration volume of sodium thiosulfate.


## 2. Apparatus

Main unit	Automatic potentiometric titrator (preamplifier STD)
Electrode	Combined platinum electrode

## 3. Reagents

Titrant	0.1mol/L sodium thiosulfate
Solvent	Pure water
Additive	10% potassium iodide

## 4. Example



—Measurement results—

	Sample (mL)	Titer (mL)	Copper (g/L)
1	10.0	18.2338	11.60
2	10.0	18.1745	11.56
3	10.0	18.2149	11.59
4	10.0	18.2019	11.58
5	10.0	18.2296	11.60
Average			11.59
SD			0.02
RSD(%)			0.1

Please feel free to contact us for any further information.  
 <Contact> Kyoto Electronics Manufacturing Co., Ltd.  
 Overseas Sales & Marketing Sect.  
<http://www.kyoto-kem.com/en/contact/form.php>