

Application Memo

Acid Number of Phenol Resin

Industry	Organic chemical industry
Instrument	Automatic potentiometric titrator
Measurement method	Acid-base titration
Standards	

1. Overview

The acid number of phenol resin is measured by titration with 0.1mol/L potassium hydroxide (ethanol) solution after the sample is dissolved in a solvent. The endpoint is determined by the maximum inflexion point on the titration curve. The acid number is calculated from the titration volume of potassium hydroxide.


2. Apparatus

Main unit	Automatic potentiometric titrator (preamplifier PTA)
Electrode	Photometric sensor Interference filter (530nm)

3. Reagents

Titrant	0.5 mol/L potassium hydroxide (ethanol)
Solvent	Ethanol, Toluene
Reagent	0.1mol/L potassium hydroxide solution
Indicator	Phenolphthalein indicator

4. Example



—Titration curve—

—Measurement results—

	Sample (g)	Titer (mL)	Acid number (mg/g)
1	2.0437	2.1654	29.839
2	2.0323	2.1670	30.029
3	2.0488	2.1797	29.962
Average			29.943
SD			0.096
RSD(%)			0.32

Please feel free to contact us for any further information.
 <Contact> Kyoto Electronics Manufacturing Co., Ltd.
 Overseas Sales & Marketing Sect.
<http://www.kyoto-kem.com/en/contact/form.php>