

Application Memo

Acidity of Tobanjan (Chinese Chili Bean Sauce)

Industry	Food & beverage
Instrument	Automatic potentiometric titrator
Measurement method	Acid base titration
Standards	

1. Overview

Acetic acid concentration as acidity of the diluted sample is measured by titration with 0.1mol/L sodium hydroxide up to pH8.3. The acetic acid concentration is calculated from the titration volume of the sodium hydroxide.


2. Apparatus

Main unit	Automatic potentiometric titrator (preamplifier STD)
Electrode	Combined glass electrode Temperature compensation electrode

3. Reagents

Titrant	0.1mol/L sodium hydroxide
---------	---------------------------

4. Example



-Titration curve-

-Measurement results-			
	Sample (g)	Titer (mL)	Concentration (mg/g)
1	2.5937	0.9436	2.2025
2	2.5432	0.9262	2.2048
3	2.5886	0.9484	2.2180
Average			2.2084
SD			0.0084
RSD (%)			0.37

Please feel free to contact us for any further information.

< Contact > Kyoto Electronics Manufacturing Co., Ltd.

Overseas Sales & Marketing Sect.

<http://www.kyoto-kem.com/en/contact/form.php>