

Application Memo

Determination of Salt in Chicken Stock Powder

Industry	Food & beverage
Instrument	Automatic potentiometric titrator
Measurement method	Precipitation titration
Standards	

1. Overview

Salt is measured by titration with 0.1mol/L silver nitrate solution. The endpoint is the maximum inflexion on the titration curve. The salt concentration is calculated from the titration volume of the silver nitrate solution.


2. Apparatus

Main unit	Automatic potentiometric titrator (preamplifier STD)
Electrode	Combined silver electrode (Internal solution: 1mol/L KNO ₃)

3. Reagents

Titrant	0.1mol/L silver nitrate solution
---------	----------------------------------

4. Example



Please feel free to contact us for any further information.
 <Contact> Kyoto Electronics Manufacturing Co., Ltd.
 Overseas Sales & Marketing Sect.
<http://www.kyoto-kem.com/en/contact/form.php>