



Application Memo Acidity of Red Wine

Industry Food & beverage

Instrument Automatic potentiometric titrator

Measurement method Acid-base titration

Standards Official Analysis Method of National Tax Agency JAPAN

1. Overview

Wine acidity is measured by neutralization titration with sodium hydroxide according to the above standard. The diluted red wine sample is titrated with 0.1ml/L sodium hydroxide up to pH8.2 and the titration volume is converted to the tartaric acid concentration as acidity. If a sample wine contains carbon dioxide gas, first boil to degas.

2. Apparatus

Main unit Automatic potentiometric titrator (preamplifier STD)


Electrode Combined glass electrode

Temperature compensation electrode

3. Reagents

Titrant 0.1mol/L sodium hydroxide

4. Example



—Measurement results—			
	Sample	Titer	Acidity
	(mL)	(mL)	(g/100mL)
1	10.00	7.0950	0.5400
2	10.00	7.0950	0.5400
3	10.00	7.0861	0.5393
Average			0.5398
SD			0.0004
RSD(%)			0.07

Please feel free to contact us for any further information.

< Contact > Kyoto Electronics Manufacturing Co., Ltd. Overseas Sales & Marketing Sect.

http://www.kyoto-kem.com/en/contact/form.php