

Application Memo

Quantitative Determination of Styrene in Latex

Industry	Plastics & rubber
Instrument	Automatic potentiometric titrator
Measurement method	Oxidation-reduction titration
Standards	

1. Overview

After adding Methanol and 15% sulfuric acid to the prepared sample, it is titrated with 0.1mol/L potassium bromate/potassium bromide solution until the color changes to yellow. The blank value is calculated from the titration volume of the titrant. After adding 10% potassium iodide, it is titrated with 0.1mol/L sodium thiosulfate solution. The endpoint is the maximum inflexion on the titration curve. The styrene concentration is calculated from the blank value and the titration volume of the sodium thiosulfate solution.


2. Apparatus

Main unit	Automatic potentiometric titrator (preamplifier STD)
Electrode	Platinum electrode Ceramic reference electrode

3. Reagents

Titrant	0.1mol/L potassium bromate/potassium bromide solution 0.1mol/L sodium thiosulfate solution
Solvent	Methanol
Additive	15% sulfuric acid, 10% potassium iodide

4. Example



—Titration curve—

—Measurement results—

	Sample size (g)	Titer KBr-KBrO ₄ (mL)	Blank value	Titer Na ₂ S ₂ O ₂ (mL)	Styrene concentration (%)
1	22.119	18.387	1.8028	0.2512	0.4179
2	22.168	18.761	1.8395	0.1647	0.4276
3	22.188	18.811	1.8444	0.1979	0.4276
4	22.189	18.753	1.8472	0.1180	0.4301
5	22.168	18.624	1.8261	0.1299	0.4253
Average					0.4257
SD					0.0047
RSD(%)					1.1

Please feel free to contact us for any further information.

<Contact> Kyoto Electronics Manufacturing Co., Ltd.

Overseas Sales & Marketing Sect.

<http://www.kyoto-kem.com/en/contact/form.php>